



MTF-Plus (Engine)

X1 17/.04

X2 25/.06

X3 30/.07

Material : Nickel-Titanium Length : 21mm,25mm,31mm

Recommended rotary speed : 250~350rpm

Package : 3pcs in one pack



Performance

Based on MTF(Engine), but MTF-Plus(Engine) with better performance and more continues after added refinements to improve performance. Rectangular cross-section design for greater strength that further enhances MTF(Engine use) canal shaping efficiency . M-Wire NiTi alloy for increased flexibility and resistance to cyclic fatigue.



Using instruction

Shaping Files (X1): Preparation the upper and middle 1/3 root canal with X1.

Finishing Files (X2, X3) All used for the formation of 1/3 root canal tip.

Application sequences

1. Establishing linear channel.
2. Using #10, #15 K file to detect the root canal.
3. X1, if necessary, V0 reaches the length of #15 K file.
4. Using #15 K file to measure and determine the working length.
5. X2 to the working length.
6. X3 to the working length.

Remark: Suggest dentist to use small Torsion force at the beginning and then can increase step by step.